## **Histocompatibility Lab | Non-Transplant Testing**

Phone: 800-245-3117 x6250 | Fax 414-937-6322



NOTE: Versiti does NOT bill patients or insurance. Test orders must be placed through a medical facility that has an account with Versiti. Client # required.

Ordering Institution Information										
Person Completing Requisition:				Physician/Provider:						
Institution:								Client #:		
Dept:			Addre	ess:						
City:				State:			Zip Code:			
Phone (Lab): Prov			der Contact (phone/email):							
Special Reporting Requests:								PO #:		
Is testing for outpatient Medicare enro		•			☐ No					
If yes, please complete the beneficiary for	rm located at https://versiti.o	org/produc	ts-serv	ices/requisition	ns and	submit wi	th this requisiti	on.		
Patient Information		F: . N						DOD		
Last Name:		First N	lame:				MI:	DOB:		
MR#:	Accessi	on #:				SSN:				
Biologic Sex/Sex Assigned at Birth:  ☐ Male ☐ Female ☐ Other	•	y: □ Ashkenazi Jewish □ Black/African American □ Central Asian □ East Asian □ Hispanic/Latino □ Middle Eastern □ Native American □ South Asian □ White □ Other								
Specimen Information										
Specimen Type: ☐ Blood ☐ Buccal Swabs ☐ Plasma ☐ Serum ☐ DNA ☐ Umbilical Cord Blood ☐ Other										
_	□ ACDB □ Clot □ Oth	ner			Draw	/ Date:		Draw 1	Γime:	
Patient History										
Potential Platelet Recipient?  Yes No If yes, please complete information below **Determine if STAT testing is required and alert Histocompatibility Lab if STAT testing is required (refer to bottom left section and drawing instructions on back/next page).										
Transfusion History:  Unknown None Multiple Last Transfusion:						/ of:				
Diagnosis: Previous Typing, If Known:										
HLA-A HLA-B HLA-C HLA-DR HLA-DQ HLA-DP										
Testing										
HLA Matched Platelet Transfusion Wo	rkup									
☐ HLA-AB Low Resolution (2303) ☐ STAT Testing (STAT Fee Applies) – complete bottom left section										
☐ HLA Antibody Identification Class I High Resolution (2226)										
HLA-Typing  ☐ HLA-DRB1 and -DQA1/-DQB: ☐ HLA-DRB1 and -DQA1/-DQB:			R1 Low	Low Resolution (2553) ☐ HLA-ARC F				C High Resolution (2329)		
☐ HLA-B Low Resolution (2305)	<u> </u>	☐ HLA-DQA1/-DQB1 Low Resolution			, ,				esolution by DNA Sequencing (2322)	
☐ HLA-C Low Resolution (2306) ☐ HLA-DPB1/-DPA1 Low Resolu				, ,		RB3*01 (DR52) Determination (5252)				
☐ HLA-ABC Low Resolution (2302) ☐ HLA-A High Resolution by DNA Sequence										
☐ HLA-DRB1 Low Resolution (2307) ☐ HLA-B High Resolution by DNA				Sequencing (2325)						
☐ HLA-DRB3, B4, B5 Low Resolution (2122) ☐ HLA-C High Resolution by DNA Sequencing (2326)										
HLA Disease Association Testing										
☐ HLA-A*02:01 Determination (2279) ☐ HLA-B51 Determination (2275)			☐ HLA-DR Single Antigen (2361) Specify:							
☐ HLA-A29 Determination (2274)	☐ HLA-B*15:02 De	☐ HLA-B*15:02 Determination (2			276) □ HLA			Typing for Narcolepsy (2270)		
☐ HLA-B27 Determination (2271)	☐ HLA-B*57:01 De	☐ HLA-B*57:01 Determination (2272)				☐ HLA Typing for Celiac Disease (2277)				
HLA Antibody Testing										
☐ HLA Antibody Detection (Flow Cytometry) (2235)  ☐ HLA Antibody Identification Class I High Resolution (2226)  ☐ HLA Antibody Identification Class II High Resolution (2231)										
STAT Testing									USE ONLY	
☐ STAT Testing (STAT Fee Applies)							HEPB	ACDA	ACDBEDTA	
Results Required No Later Than: Da	te Needed By:/		Time	:::_			Clot	Other: _		
Contact Name:							Opened By:		Reviewed By:	
Contact Phone #:							Evaluated By:		Labeled By:	

NOTE: Versiti does NOT bill patients or insurance. Test orders must be placed through a medical facility that has an account with Versiti. Client # required.

## **DRAWING INSTRUCTIONS**

Tubes must be <u>individually</u> labeled with **FULL NAME OF INDIVIDUAL, ANOTHER IDENTIFIER** (e.g., SSN, MRN, DOB), DATE AND TIME OF DRAWING. Samples cannot be accepted after any exposure to an environment in which HLA genes are amplified. This precaution is essential to avoid contamination of samples with DNA that could alter test results. Samples will be accepted from 8:00 a.m. Monday through noon on Friday. Emergency testing MUST be arranged through the laboratory. Call (414) 937-6201.

TEST	SAMPLE REQUIREMENTS	STORE and SHIP
HLA Low Resolution (A, B, C, AB, ABC, DRB1, DRB3 B4 B5, DQB1/DQA1, DPB1), HLA-A29, HLA-B27, HLA-B51, HLA-Narcolepsy, HLA-Celiac, HLA-DR Single Antigen	5 - 14 ml EDTA whole blood (lavender top) <b>or</b> ACDA whole blood (yellow top) or 4 buccal swabs	Room temperature
HLA High Resolution (A, B, C, DRB1, DQB1, DPB1), HLA-A*02:01, HLA-B*15:02, HLA-B*57:01	5 - 14 ml EDTA whole blood (lavender top) <b>or</b> ACDA whole blood (yellow top) or 4 buccal swabs	Room temperature
HLA Antibody Detection, HLA Antibody Identification	10 ml Clotted (red top) blood	Room temperature
Platelet Recipient: HLA-AB LowResolution HLA Antibody Identification Class I High Resolution * Indicate STAT Testing, if required	14 ml EDTA whole blood (lavender top) <b>or</b> 4 buccal swabs <b>and</b> 10 ml Clotted (red top) blood	Room temperature
Kidney, Heart, Liver, Pancreas, Lung Recipient - Initial Workup	REFER TO HISTOCOMPATIBILITY LAB TRANSPLANT TESTING REQUISITION	Room temperature
Kidney Donor Workup	REFER TO HISTOCOMPATIBILITY LAB TRANSPLANT TESTING REQUISITION	Room temperature
Crossmatch (Flow Cytometric)	REFER TO HISTOCOMPATIBILITY LAB TRANSPLANT TESTING REQUISITION	Room temperature
Bone Marrow (Stem Cell) recipients or donors	REFER TO HISTOCOMPATIBILITY LAB TRANSPLANT TESTING REQUISITION	Room temperature

## **SHIPPING INFORMATION**

Contact laboratory for pediatric drawing requirements or low white cell count drawing requirements. Blood samples should be shipped by overnight carrier. The package must be shipped in compliance with carrier's guidelines. Please contact your carrier for current biohazardous shipping regulations.

Label Box: Refrigerate, Room Temperature, or Frozen (whichever is appropriate)

Packages should be addressed to:

Versiti Wisconsin – Histocompatibility Laboratory 638 N 18<sup>th</sup> Street Milwaukee, WI 53233